



# OXLIFE INDEPENDENCE<sup>®</sup>

## YOUR OXYGEN PRESCRIPTION AT HOME OR OUTDOORS

- 12 pulse settings and 3 LPM continuous flow, giving you total flexibility to comfortably travel \*Best in Class\*
- Swappable batteries allowing you to extend the battery time without needing to charge the entire device
- 2 in 1: Functions as a both a portable and stationary device in only 1 device
- Robust, reliable, and durable with a sturdy frame and inbuilt wheels
- Transportable: Easy to carry in a car, train, or plane
- Longest battery run time \*Best in Class\*



2 year warranty



Made in the US



FAA approved for  
Air Travel



## SUPERIOR PRODUCT FEATURES

- Industry Leading Durability and Reliability = Lower Operating Costs
- Energy Smart Technology® = Longest Battery Life in Class
- Easy Battery Swap While Device is in Operation
- Lightweight Metal / Polycarbonate with **\*Best in Class\*** Components
- Integrated Cart with 6-inch Wheels Offers **\*Best in Class\*** Maneuverability

## PRODUCT SPECIFICATIONS

Dimensions	Height: 20.2" • Width: 10.9" • Depth: 9.45"
Weight	7.5 kgs (with cart & battery)
Battery Recharge Time	2.5 hours
Oxygen Concentration	87% - 95% at sea level at all settings
DC Power	12 - 15 VDC • 15 amps
AC Power	Universal Input: 100 - 240 VAC 50 - 60 Hz • 2.0A maximum
Battery Duration	5.75 hours @ 2 pulse at 20 bp 2.5 hours @ 2 LPM continuous
Oxygen Flow Rate (Pulse 15-40 BPM)	Continuous: 0.5 - 3 LPM in 0.5 increments Pulse: 0.5 - 6.0 in 0.5 increments (8 - 96 ML)

**Customisable:** Additional batteries, car charger, and external battery charger are available for extended travel



### Marketed by

Office No. 506, 5th Floor, Office Tower – B, Elante,  
Plot no. 178-178A, Ind Area, Phase I, Chandigarh

Call : +91 7341103371  
[www.sanraishop.com](http://www.sanraishop.com)