


OxyPure

- 1-5L/min
- Fast & Easy Maintenance
- Lightweight and Robust Design
- Oxygen and Performance Monitor
- Suitable for Most Oxygen Prescriptions



OxyPure

OXYGEN CONCENTRATOR



Specification	Model: Oxypure
Oxygen Flow	0.5 to 5 Liters per minute
Delivery Method	Continuous
Sound Levels (dBA)	≤40
Weight (kg)	15.2
Size (LxWxH)	370 x 330 x 580mm
Cord Length (m)	2 m
Oxygen Purity	93% ± 3%
Oxygen Outlet Pressure	8.5psi ± 0.5psi
Electrical Requirements	220V/50Hz



Marketed by

Office No. 506, 5th Floor, Office Tower – B, Elante,
Plot no. 178-178A, Ind Area, Phase I, Chandigarh

Call : +91 7341103371
www.sanraishop.com